



# Efficient cleanup of low molecular weight carbohydrate samples using gravity flow columns

**Intellectual Property Notice:** The Biopharma business of GE Healthcare was acquired by Danaher on 31 March 2020 and now operates under the Cytiva™ brand. Certain collateral materials (such as application notes, scientific posters, and white papers) were created prior to the Danaher acquisition and contain various GE owned trademarks and font designs. In order to maintain the familiarity of those materials for long-serving customers and to preserve the integrity of those scientific documents, those GE owned trademarks and font designs remain in place, it being specifically acknowledged by Danaher and the Cytiva business that GE owns such GE trademarks and font designs.

## cytiva.com

GE and the GE Monogram are trademarks of General Electric Company. Other trademarks listed as being owned by General Electric Company contained in materials that pre-date the Danaher acquisition and relate to products within Cytiva's portfolio are now trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva. Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. All other third-party trademarks are the property of their respective owners.  
© 2020 Cytiva  
All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.  
For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

# Efficient Cleanup of Low Molecular Weight Carbohydrate Samples using Gravity Flow Columns

Ola Rönn <sup>1\*</sup>, Helena Hedlund <sup>1</sup>, and Dorothe Spillmann <sup>2</sup>

<sup>1</sup>GE Healthcare Bio-Sciences AB, Björkgatan 30, SE-751 84 Uppsala, Sweden. \*ola.ronn@ge.com

<sup>2</sup>Dept. of Medical Biochemistry and Microbiology, Uppsala University, SE-751 23 Uppsala, Sweden.

## Introduction

Low molecular weight carbohydrates (oligosaccharides) can be difficult to separate from contaminants such as salts. Methods like ion exchange are not suitable due to interference by the contaminant. Sephadex™ G-10 is a well established gel filtration chromatography medium suitable for group separation of biomolecules with a molecular weight above 700 from smaller molecules such as salts, dyes, and free radioactive isotopes. Here, we show efficient cleanup of an oligosaccharide-containing solution from salt using PD MiniTrap™ G-10, a new pre-packed column format.



## Conclusions

- Rapid cleanup of carbohydrates with molecular weights above 700 prior to downstream analysis
- Oligosaccharide recovery of ~95%
- Efficient removal of contaminants such as salts, dyes, and radioactive labels
- High desalting capacity



## Methods

### Sample preparation

- $^3\text{H}$ -labeled octasaccharides were prepared from heparin by deaminative cleavage, reduced with  $\text{NaB}_3\text{H}_4$  followed by fractionation on a size exclusion column.
- Concentrated  $^3\text{H}$ -labeled oligo- or polysaccharides in water were diluted with different salt solutions to the final concentrations of salt as indicated.

### Salt removal using PD MiniTrap G-10 gravity column

- The column was prepared for use by resuspending the medium, letting it settle, and then the storage solution was removed by gravity flow.
- The column was washed with five column volumes of distilled water.
- A sample of 100  $\mu\text{l}$   $^3\text{H}$ -heparin octamer (20 000 cpm) in 0.5 M NaCl was applied on top of the gel bed and allowed to move into the gel. (elution profile start, Fig. 1)
- After sample entry the column was eluted with 20x100  $\mu\text{l}$ -aliquots of water with fraction collection of equal aliquots.

### Analysis of eluate


- One milliliter of Optiphase™ HiSafe® 3 scintillation cocktail was added to an aliquot of each sample (50  $\mu\text{l}$ ) and analyzed in a Beckman Coulter® LS6500 multipurpose scintillation counter.
- Conductivity was measured by diluting 50  $\mu\text{l}$  of eluate to 1000  $\mu\text{l}$  with distilled water and analyzing on a Radiometer CDM 83 conductivity meter.

## Results

Efficient cleanup of the carbohydrate solution was achieved with recovery and salt removal typically over 90%. The elution profile for the octasaccharide showed excellent peak separation from the salt contaminant (Fig. 1).

Corresponding experiments using peptides as the target gave similar results (not shown). Analyses also show that the limits of separation is set by sample volume rather than concentration. For sample volumes exceeding 300  $\mu\text{l}$ , a larger column format is needed, for example PD MidiTrap™ G-10.

### Elution profile



**Fig 1.** Removal of NaCl from a  $^3\text{H}$ -heparin octasaccharide solution on a PD MiniTrap G-10 column. The fractions were analyzed with regard to sample content (radioactivity; blue) and salt content (conductivity, red). The oligosaccharide recovery was approximately 95% and the desalting capacity was 94% (between arrows).

## Sephadex G-10 characteristics

Matrix	Cross-linked dextran
Separation mechanism	According to size
Particle size range	55-165 $\mu\text{m}$
Exclusion limit ( $M_r$ )	700
Chemical stability	All common buffers
Working pH range	2-13

### Characteristics of prepacked columns

	MiniTrap G-10	MidiTrap G-10
Vol. prepacked medium	2.1 ml	5.3 ml
Packed bed dimensions	0.97 x 2.8 cm	1.3 x 4.0 cm
Column volume	5 ml	8.5 ml
Void volume	700 $\mu\text{l}$	1.6 ml
Max. sample volume	300 $\mu\text{l}$	1.0 ml
Vol. eluted sample	0.5 ml	1.2 ml
Recovery <sup>1</sup>	70% to 90%	70% to 90%
Desalting capacity	>75%	>75%
Column material	PP and PE	PP and PE
Storage solution	20% ethanol	20% ethanol
Storage temp.	4°C to 30°C	4°C to 30°C

<sup>1</sup> Biomolecule dependent

**Fig 2.** PD MiniTrap G-10 columns are prepacked with Sephadex G-10. The maximum sample volume is 300  $\mu\text{l}$ . For sample volumes up to 1 ml PD MidiTrap G-10 columns are also available.

## Ordering information

Product	Quantity	Code No.
PD MiniTrap G-10	50 columns	28-9180-10
PD MidiTrap G-10	50 columns	28-9180-11

Literature	Code No.
Gel Filtration Principles and Methods, Handbook	18-1022-18
PD MiniTrap G-10, PD MidiTrap G-10 data file	28-9267-49

For contact information for your local office,  
please visit, [www.gelifesciences.com/contact](http://www.gelifesciences.com/contact)

[www.gelifesciences.com/protein-purification](http://www.gelifesciences.com/protein-purification)

GE Healthcare Bio-Sciences AB  
Björkgatan 30  
751 84 Uppsala  
Sweden

GE, imagination at work and GE monogram are trademarks of General Electric Company.  
Sephadex, MiniTrap, and MidiTrap are trademarks of GE Healthcare companies.

All third party trademarks are the property of their respective owners.

© 2008 General Electric Company – All rights reserved.

First published Apr. 2008

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

