



Application Note

USD 3051

High Viscosity Mixing with a Pall® Wand Mixer

Mixing system: Pall Wand Mixer

Mixing biocontainer: 200 L Mixer biocontainer

Application mixing type: Liquid-liquid

The Pall Wand Mixer is a compact and non-invasive single-use mixing system. The heart of this system is a mixing biocontainer incorporating an innovative top-mounted impeller capable of providing efficient mixing for powder-liquid and liquid-liquid mixing applications. The impeller comprises a rotating wand inside an inert polymer sleeve, and is designed to ensure low particle shedding and total containment while serving effectively in a wide variety of mixing tasks.

Introduction

High-viscosity materials can present a formidable mixing challenge, yet they also serve to illustrate the relative strengths and weakness of different mixing technologies. In this experiment, a Wand Mixer was used to perform liquid-liquid mixing in a highly viscous matrix – high fructose corn syrup 90% at a viscosity of 1200 centipoise.

Experimental

A 200 L Wand Mixer mixing biocontainer was filled with 200 L of ADM Cornsweet♦ 90 high fructose corn syrup (HFCS90), and mixing speed was set to 250 rpm. A 60 mL aliquot of dilute sodium hydroxide solution (8 M NaOH) was then added, via a top-mounted inoculation port, to the surface of the HFCS90. Solution homogeneity was monitored via a pH sensor mounted in the top of the biocontainer.




Results

Figure 1 shows solution homogeneity in the biocontainer during mixing. After the NaOH addition, mixing was found to be complete within approximately 5 minutes.

Figure 1

Shows solution homogeneity in the biocontainer during mixing



Conclusions

The Pall Wand Mixer can perform effectively at viscosity up to 1200 cP. For mixing applications <200 L maximum volume, the 5 minute mixing time in this experiment makes the Pall Wand Mixer an excellent choice when high-viscosity liquid-liquid mixing is the primary intended purpose.



Life Sciences

Corporate Headquarters

Port Washington, NY, USA
+1.800.717.7255 toll free (USA)
+1.516.484.5400 phone
biopharm@pall.com e-mail

European Headquarters

Fribourg, Switzerland
+41 (0)26 350 53 00 phone
LifeSciences.EU@pall.com e-mail

Asia-Pacific Headquarters

Singapore
+65 6389 6500 phone
sgcustomerservice@pall.com e-mail


Visit us on the Web at www.pall.com/mixing

E-mail us at allegro@pall.com

International Offices

Pall Corporation has offices and plants throughout the world in locations such as: Argentina, Australia, Austria, Belgium, Brazil, Canada, China, France, Germany, India, Indonesia, Ireland, Italy, Japan, Korea, Malaysia, Mexico, the Netherlands, New Zealand, Norway, Poland, Puerto Rico, Russia, Singapore, South Africa, Spain, Sweden, Switzerland, Taiwan, Thailand, the United Kingdom, the United States, and Venezuela. Distributors in all major industrial areas of the world. To locate the Pall office or distributor nearest you, visit www.pall.com/contact.

The information provided in this literature was reviewed for accuracy at the time of publication. Product data may be subject to change without notice. For current information consult your local Pall distributor or contact Pall directly.

© 2015, Pall Corporation. Pall and  are trademarks of Pall Corporation. ® indicates a trademark registered in the USA and TM indicates a common law trademark. *Filtration.Separation.Solution.* is a service mark of Pall Corporation. ♦Cornsweet is a trademark of Archer Daniels Midland Company.