



Convenient purification of monoclonal antibodies using HiTrap rProtein A FF

Intellectual Property Notice: The Biopharma business of GE Healthcare was acquired by Danaher on 31 March 2020 and now operates under the Cytiva™ brand. Certain collateral materials (such as application notes, scientific posters, and white papers) were created prior to the Danaher acquisition and contain various GE owned trademarks and font designs. In order to maintain the familiarity of those materials for long-serving customers and to preserve the integrity of those scientific documents, those GE owned trademarks and font designs remain in place, it being specifically acknowledged by Danaher and the Cytiva business that GE owns such GE trademarks and font designs.

cytiva.com

GE and the GE Monogram are trademarks of General Electric Company.
Other trademarks listed as being owned by General Electric Company contained in materials that pre-date the Danaher acquisition and relate to products within Cytiva's portfolio are now trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.
Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.
All other third-party trademarks are the property of their respective owners.
© 2020 Cytiva
All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.
For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

Convenient Purification of Monoclonal Antibodies using HiTrap rProtein A FF

M. Carlsson, A. Heijbel, and A-C. Häggqvist

GE Healthcare AB, Björkgatan 30, SE-751 84 Uppsala, Sweden

Abstract

Purification of monoclonal antibodies can be achieved by a number of methods, protein A affinity chromatography being especially powerful. HiTrap™ rProtein A FF is a prepacked column for the purification of monoclonal antibodies in a quick, easy, and convenient way. The work presented shows different applications using HiTrap rProtein A FF 1 ml and 5 ml columns. Mouse monoclonal antibodies were purified from cell cultures and ascitic fluids. The purifications were performed using either a syringe or chromatography system.

Humanized IgG₄ was purified directly from a myeloma cell culture resulting in a highly purified antibody at a yield of 93%. The capacity for a mouse IgG_{2a} purified from cell culture supernatant was 18 mg/ml medium at a flow rate of 624 cm/h (4 ml/min on a 1 ml column), showing the high capacity of HiTrap rProtein A FF at high flow rates. All purified monoclonals were >95% pure according to analysis by SDS-PAGE and silver staining. The optimum binding conditions for a mouse IgG₁ were investigated in a series of scouting experiments performed on ÄKTAexplorer™. This particular antibody was efficiently bound to HiTrap rProtein A FF in a phosphate buffer containing 0.15 M NaCl at neutral pH.

Introduction

Purification of monoclonal antibodies is required for many applications in numerous fields of science and technology. Protein A affinity chromatography is well known as one of the most efficient techniques for the purification of monoclonal antibodies. Due to the high specificity of protein A for immunoglobulins, high product purity is achieved within one chromatographic step.

A rapid and convenient way to purify monoclonal antibodies is to use HiTrap rProtein FF columns. This prepacked, ready to use column, 1 ml or 5 ml, contains rProtein A Sepharose™ Fast Flow. The ligand, recombinant protein A, has been specially engineered with a C-terminal cysteine to allow the formation of a stable thioether link between the ligand and the base matrix. Terminal coupling allows the ligand to exert maximum binding capacity for IgG.

Using HiTrap rProtein A FF, monoclonal antibodies from cell cultures and ascitic fluids can easily be prepared to high purity in a single step. The column can either be operated by a simple syringe, or by coupling it to a chromatographic system. To achieve even higher binding capacities in one run, several HiTrap rProtein A FF columns can be connected together.

Conclusions

- HiTrap rProtein A FF provides rapid and convenient preparative purification of monoclonal antibodies from cell cultures and ascitic fluids.
- High product purity (>95%) is achieved within one chromatographic step.
- High binding capacities are obtained even at high flow rates.
- HiTrap rProtein A FF can be operated with a simple syringe or by connecting it to a chromatographic system, such as ÄKTAexplorer.



Material & Methods

Columns:	HiTrap rProtein A FF, 1 ml or 5 ml
Samples:	Different cell culture supernatants or ascitic fluids containing different subclasses of IgG
Sample pre-treatment:	Filtered through a Millipore Millex® 0.45 µm filter
Binding buffer:	0.02 M sodium phosphate, pH 7.0 (incl. 3 M NaCl for mouse IgG ₁)
Elution buffer:	0.1 M sodium citrate, pH 3.0
Neutralizing buffer:	1 M Tris-HCl, pH 9.0
Fraction treatment:	200 µl neutralizing buffer/ml collected eluted fraction

Characteristics of HiTrap rProtein A FF

Column dimensions, i.d. x h:	0.7 x 2.5 cm (1 ml) 1.6 x 2.5 cm (5 ml)
Ligand:	recombinant protein A, (<i>E. coli</i>)
Coupling chemistry:	epoxy
Degree of substitution:	≈ 6 mg rProtein A/ml medium
Total binding capacity:	≈ 50 mg human IgG/ml medium
Dynamic binding capacity for some monoclonal antibodies*:	23 mg mouse IgG _{2b} /ml medium 12 mg mouse IgG ₁ /ml medium 11 mg humanised IgG ₄ /ml medium
Average particle size:	90 µm
Bead structure:	highly cross-linked 4% agarose
Maximum backpressure:	3 bar, 0.3 MPa
Maximum flow rate at room temperature:	4 ml/min (624 cm/h) and 20 ml/min (600 cm/h) for 1 and 5 ml columns, respectively
Recommended flow rate:	1 ml/min (156 cm/h) and 5 ml/min (150 cm/h) for 1 and 5 ml columns, respectively
Chemical stability:	All commonly used buffers
pH stability:	
Working	3–10
Cleaning	2**–11
Temperature stability:	
Working	4°C–40°C
Storage	4°C–8°C
Storage buffer:	20% ethanol

* Running conditions:


Column:	HiTrap rProtein A FF 1 ml
Sample:	cell culture supernatants
Binding buffer:	20 mM sodium phosphate (incl. 3 M NaCl for IgG1), pH 7.0.
Elution buffer:	0.1 M sodium citrate, pH 3.0
Flow rate:	1 ml/min.

** pH below 3 is sometimes required to elute strongly bound Ig's. However, protein ligands may hydrolyze at very low pH.

Applications


Purification of humanized IgG₄ from a myeloma cell culture

Column:	HiTrap rProtein A FF 1 ml
Sample volume:	30 ml containing 12 mg IgG ₄ The sample was a kind gift from Dr. J. Bonnerjea, LONZA Biologics plc, UK
Binding buffer:	0.02 M sodium phosphate, pH 7.0
Elution buffer:	0.1 M sodium citrate, pH 3.0
Flow rate:	1 ml/min (156 cm/h)



SDS-PAGE

Lane 1:	Low Molecular Markers
Lane 2:	Starting material
Lane 3:	Flow-through pool
Lane 4:	Eluted IgG ₄ pool



Result


Eluted amount of IgG₄: 11.2 mg

Yield: 93%

Purity: > 95% according to SDS-PAGE

Purification of mouse IgG₁ from ascitic fluid using a syringe

Column: HiTrap rProtein A FF 1 ml
 Sample volume: 1 ml (0.5 ml ascites diluted with 0.5 ml binding buffer)
 The sample was a kind gift from Dr. N. Linde, EC Diagnostics, Sweden
 Binding buffer: 0.02 M sodium phosphate, 3 M NaCl, pH 7.0
 Elution buffer: 0.1 M sodium citrate, pH 3.0
 Flow rate: ≈ 1 ml/min
 Instrumentation: Syringe




SDS-PAGE

Lane 1: Low Molecular Markers
 Lane 2: Starting material, dil. 20x
 Lane 3: Flow-through pool, dil. 5x
 Lane 4: Eluted IgG₁ pool, dil. 5x
 Lane 5: Flow-through pool, dil. 5x, reduced
 Lane 6: Eluted IgG₁ pool, reduced
 Lane 7: Low Molecular Markers


Purification of mouse IgG_{2a} from cell culture supernatant

Column: HiTrap rProtein A FF 1 ml
 Sample volume: 150 ml
 Binding buffer: 0.02 M sodium phosphate, pH 7.0
 Elution buffer: 0.1 M sodium citrate, pH 3.0
 Flow rate: 4 ml/min (624 cm/h)




SDS-PAGE

Lane 1: Low Molecular Markers reduced
 Lane 2: Starting material, reduced
 Lane 3: Eluted IgG_{2a} pool, dil. 3x, reduced



Result

Eluted amount of IgG₁: 1.1 mg
 Purity: > 95% according to SDS-PAGE




Result

Eluted amount of IgG_{2a}: 18.2 mg
 Purity: > 95% according to SDS-PAGE


Purification of mouse IgG_{2a} from ascitic fluid

Column: HiTrap rProtein A FF 1 ml
Sample volume: 7 ml (3.5 ml ascites diluted with 3.5 ml binding buffer)
Binding buffer: 0.02 M sodium phosphate, 3 M NaCl, pH 7.0
Elution buffer: 0.1 M sodium citrate, pH 3.0
Flow rates:
 sample application 0.5 ml/min (78 cm/h)
 wash and elution 1 ml/min (156 cm/h)



SDS-PAGE

- Lane 1:** Low Molecular Markers, reduced
Lane 2: Starting material, dil. 40X, reduced
Lane 3: Starting material, dil. 20X, reduced
Lane 4: Flow-through pool, dil. 20X, reduced
Lane 5: Flow-through pool, dil. 10X, reduced
Lane 6: Eluted IgG_{2a} pool, dil. 10X, reduced
Lane 7: Eluted IgG_{2a} pool, dil. 5X, reduced
Lane 8: Low Molecular Markers, reduced



Result

Eluted amount of IgG_{2a}: 14.5 mg
 Purity: > 95% according to SDS-PAGE

Scouting for binding conditions using ÄKTAexplorer with BufferPrep*



The concentration of NaCl in the binding buffer was varied over the range 0.15 M to 2.5 M at constant pH.

Column: HiTrap rProtein A FF 5 ml
Sample: Cell culture supernatant containing mouse IgG₁
Sample volume: 8 ml
Flow rate: 2 ml/min (60 cm/h)
Instrumentation: ÄKTAexplorer

BufferPrep solutions:
 1. *Buffer stock* 0.1 M sodium phosphate, 0.1 M sodium citrate
 2. *Acid stock* 0.2 M HCl
 3. *Salt stock* 5 M NaCl
 4. *Water*


Salt concentrations in the binding buffers: 0.15, 0.5, 1.0, 1.5, 2.0 and 2.5 M. pH of binding buffer 7.4

* A. Holtz, Phadia, Sweden is gratefully acknowledged for performing these experiments.



SDS-PAGE and Western Blotting from the 0.15 M NaCl run

Lane 1: Flow-through fraction 2
Lane 2: Flow-through fraction 2, reduced
Lane 3: Eluate, fraction 9
Lane 4: Eluate, fraction 10
Lane 5: Eluate, fraction 10, reduced
Lane 6: Low Molecular Markers



Result

The same amount of this particular IgG₁ bound at all NaCl concentrations tested.

SDS-PAGE analysis of the eluted fractions from the run, using 0.15 M NaCl, shows that the purified monoclonal is over 95% pure. Western blotting confirms that no antibody is in the flow-through material.

Analysis

Concentration determination

Absorbance measurement at 280 nm using Ultrospec Plus:

$$A_{280}^{1\%} = 14$$

Purity check:

SDS-PAGE electrophoresis

Gel:	PhastGel™ Gradient 10–15,
Sample pretreatment:	Dilution 1:5 with 15% SDS, 60 mM Tris, 6 mM EDTA, 0.06% Bromophenol Blue, pH 8.0 2-mercaptoethanol (final conc. 6%) was added when run under reducing conditions Heating, 5 min., 90°C
Sample volume:	1 µl
Molecular weight standard:	Low Molecular Markers
Staining:	Silver, according to the standard protocol for PhastSystem™
Instrumentation:	PhastSystem

www.gelifesciences.com/hitrap
www.gelifesciences.com/protein-purification

GE Healthcare Bio-Sciences AB
Björkgatan 30
751 84 Uppsala
Sweden

GE, imagination at work, and GE monogram are trademarks of General Electric Company.

ÅKTAexplorer, Drop Design, HiTrap, PhastGel, PhastSystem, and Sepharose are trademarks of GE Healthcare companies.

All third party trademarks are the property of their respective owners.

© 1997-2007 General Electric Company - All rights reserved.
First published 1997.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare Europe GmbH
Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare UK Ltd
Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Bio-Sciences Corp
800 Centennial Avenue, P.O. Box 1327, Piscataway,
NJ 08855-1327, USA

GE Healthcare Bio-Sciences KK
Sanken Bldg. 3-25-1, Hyakunincho, Shinjuku-ku,
Tokyo 169-0073, Japan

Asia Pacific T +85 65 62751830 F +85 65 62751829 • Australasia T +61 2 8820 8299 F +61 2 8820 8200 • Austria T 01 /57606 1613 F 01 /57606 1614 • Belgium T 0800 73 890 F 02 416 8206 • Canada T 1 800 463 5800 F 1 800 567 1008 • Central & East Europe T +43 1 972 720 F +43 1 972 722 750 • Denmark T +45 70 25 24 50 F +45 45 16 2424 • Eire T 1 800 709992 F +44 1494 542010 • Finland & Baltics T +358 9 512 3940 F +358 9 512 39439 • France T 01 69 35 67 00 F 01 69 41 98 77 Germany T 0800 9080 711 F 0800 9080 712 • Greater China T +852 2100 6300 F +852 2100 6338 • Italy T 02 26001 320 F 02 26001 399 • Japan T 81 3 5331 9336 F 81 3 5331 9370 • Korea T 82 2 6201 3700 F 82 2 6201 3803 • Latin America T +55 11 3933 7300 F +55 11 3933 7304 • Middle East & Africa T +30 210 96 00 687 F +30 210 96 00 693 • Netherlands T 0800-82 82 82 1 F 0800-82 82 82 4 • Norway T +47 815 65 777 F +47 815 65 666 • Portugal T 21 417 7035 F 21 417 3184 • Russia, CIS & NIS T +7 495 956 5177 F +7 495 956 5176 Spain T 902 11 72 65 F 935 94 49 65 • Sweden T 018 612 1900 F 018 612 1910 • Switzerland T 0848 8028 10 F 0848 8028 11 • UK T 0800 515 313 F 0800 616 927 • USA T +1 800 526 3593 F +1 877 295 8102

