



Depyrogenation of buffer or purified water

Intellectual Property Notice: The Biopharma business of GE Healthcare was acquired by Danaher on 31 March 2020 and now operates under the Cytiva™ brand. Certain collateral materials (such as application notes, scientific posters, and white papers) were created prior to the Danaher acquisition and contain various GE owned trademarks and font designs. In order to maintain the familiarity of those materials for long-serving customers and to preserve the integrity of those scientific documents, those GE owned trademarks and font designs remain in place, it being specifically acknowledged by Danaher and the Cytiva business that GE owns such GE trademarks and font designs.

cytiva.com

GE and the GE Monogram are trademarks of General Electric Company.
Other trademarks listed as being owned by General Electric Company contained in materials that pre-date the Danaher acquisition and relate to products within Cytiva's portfolio are now trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.
Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.
All other third-party trademarks are the property of their respective owners.
© 2020 Cytiva
All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.
For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

CY14305-09Jun20-AN

Depyrogenation of buffer or purified water

The situation

Pyrogens, substances that cause fever, must be removed from purified water and a host of buffers and small molecular solutions used in pharmaceutical research and production. Pyrogens take two forms, exogenous and endogenous. For the most part, exogenous pyrogens are microbes, microbial products or toxins, such as gram-negative endotoxin (bacterial lipopolysaccharides [LPS]) or toxins from *Staphylococcus* strains. Endogenous pyrogens are cytokines (polypeptides [PPs]) produced generally by monocyte macrophages but also by skin, endothelial, epithelial, and glial cells.

Maximum safety is enhanced by depyrogenating solutions immediately prior to use. Simple sanitization by heat does not remove the toxic compounds associated with LPS or PPs. Sterilization via bactericides or other chemicals results in unwanted chemical residue. A simple, consistent means of removing pyrogens at the point of use—one that does not permit last-minute contamination—is mandatory for many R&D and medical settings.

The solution: GE Healthcare hollow fiber ultrafiltration cartridges

One of the most effective ways of water depyrogenation is via membrane separations, which removes contaminants and pyrogens at the molecular level using an easily-configured system. Cross flow membrane separations lend themselves to closed systems that provide good protection against contamination.

General guidelines

- Suited to depyrogenation of purified water, buffers, and solutions of small molecules
- Pyrogens can be endotoxins or other entities
- Endotoxins (LPS) readily aggregate to a molecular weight (MW) of 250,000 kD+
- Subunits are 10 kD
- Ultrafilters can reduce endotoxin levels > 4 logs
- Sanitize and depyrogenate equipment beforehand
- Confirm with WFI rinse
- Sanitize frequently
- Maintain flow in water systems—even when not in use (Fig. 2)

Pyrogen removal from municipal water using a 10,000 nominal molecular weight cutoff (NMWC) hollow fiber ultrafiltration membrane

Membrane: 10,000 NMWC

Model number: GE Healthcare UFP-10-C-85

Feed: 0.2 µm filtered tap water

Inlet pressure: 1.03 barg (15 psig)

Outlet pressure: 0.97 barg (14 psig)

Time (minutes)	Feed	Bleed	Reject	Permeate
0	93	93	104	0.023
30	119	119	183	0.031
60	125	125	189	0.039
90	155	155	268	0.049
120	159	296	-	0.02
150	159	304	-	0.034
180	157	327	-	0.019

Figure 1. Typical depyrogenation results (in eu/ml) showing a consistently low level in the filtrate despite a 3x increase in the LPS level in the feed



Depyrogenation of USP water




Figure 2. Typical depyrogenation system utilizing ultrafiltration. This is a closed system, with no element exposed to the environment. Water from a reverse osmosis/deionization (RO/DI) system is fed to a reservoir, which then is circulated through a hollow fiber ultrafiltration cartridge. Water flows continuously, both via a reject line in the retentate circuit and a shunt in the permeate circuit that recirculates permeate when it is not being drawn off. Valves and pressure gauge are shown; necessary pumps are not.

www.gehealthcare.com

Global Headquarters GE Healthcare
Little Chalfont
Buckinghamshire, U.K. HP7 9NA

GE, imagination at work and GE monogram are trademarks of General Electric Company.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. GE Healthcare reserves the right, subject to any regulatory and contractual approval, if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation.

© 2006 General Electric Company - All rights reserved.

GE Healthcare Bio-Sciences AB, a General Electric Company.

GE Healthcare Bio-Sciences AB
Björkgatan 30, 751 84 Uppsala, Sweden

GE Healthcare Europe GmbH
Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare UK Ltd
Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Bio-Sciences Corp
800 Centennial Avenue, P.O. Box 1327
Piscataway, NJ 08855-1327, USA

GE Healthcare Bio-Sciences KK
Sanken Bldg. 3-25-1, Hyakunincho, Shinjuku-ku,
Tokyo 169-0073, Japan

Asia Pacific Tel +65 6275 1830 Fax +65 6275 1829 Australasia Tel +61 2 9899 0999 Fax +61 2 9899 7511 Austria Tel 01/57606-1619 Fax 01/57606-1627 Belgium Tel 0800 73 888 Fax 02 416 82 06 Canada Tel 800 463 5800 Fax 800 567 1008
Central, East, & South East Europe Tel +43 1 972720 Fax +43 1 97272 2750 Denmark Tel 45 16 2400 Fax 45 16 2424 Finland & Baltics Tel +358 (0)9 512 39 40 Fax +358 (0)9 512 39 439 France Tel 01 69 35 67 00 Fax 01 69 41 96 77
Germany Tel 089 96281 660 Fax 089 96281 620 Greater China Tel +852 2100 6300 Fax +852 2100 6338 Italy Tel 02 27322 1 Fax 02 27302 212 Japan Tel +81 3 5331 9336 Fax +81 3 5331 9370 Latin America Tel +55 11 3933 7304
Middle East & Africa Tel +30 210 9600 687 Fax +30 210 9600 693 Netherlands Tel 0800 82 82 82 1 Fax 0800 82 82 82 4 Norway Tel 815 65 555 Fax 815 65 666 Portugal Tel 21 417 7035 Fax 21 417 3184 Russia & other C.I.S. & N.I.S Tel +7 (495) 956 5177
Fax +7 (495) 956 5176 Spain Tel 93 594 49 50 Fax 93 594 49 55 Sweden Tel 018 612 1900 Fax 018 612 1910 Switzerland Tel 0848 8028 12 Fax 0848 8028 13 UK Tel 0800 616928 Fax 0800 616927 USA Tel 800 526 3593 Fax 877 295 8102



imagination at work