

GSTPrep FF 16/10

Instructions for use

Introduction

GSTPrep™ FF 16/10 is a prepacked, ready to use column containing Glutathione Sepharose™ 4 Fast Flow designed for one-step purification of glutathione S-transferase (GST) tagged proteins, other glutathione S-transferases, and glutathione binding proteins.

Column data

Matrix	Glutathione Sepharose 4 Fast Flow
Bead structure	Highly cross-linked 4% agarose
Mean particle size	90 µm
Ligand	Glutathione
Ligand concentration	120–320 µmol glutathione/ml medium
Binding capacity ¹	= 10 mg recombinant glutathione S-transferase/ml medium (GST, Mr 26 000)
Bed volume	20 ml
Bed height	100 mm
i.d.	16 mm
Column hardware	Polypropylene
Recommended flow rate ¹²	1–10 ml/min (30–300 cm/h)
Maximum flow rate ¹²	10 ml/min (300 cm/h)
Maximum pressure over the packed bed during operation, Δp	0.15 MPa, 1.5 bar, 22 psi
HiPrep™ column hardware pressure limit ³	0.5 MPa, 5 bar, 73 psi
pH stability	pH 3–12
Storage	4°C to 30°C in 20% ethanol

¹ Note: The binding of GST to glutathione is depending on the flow rate and therefore can often a low flow rate increase the binding capacity. This is important during loading of sample and elution.

² Water at room temperature. Flow rate is determined by $v \cdot \eta < 10 \text{ ml/min}$ where v = flow rate and η = viscosity.

³ Many chromatography systems are equipped with pressure gauges to measure the pressure at a particular point in the system, usually just after the pumps. The pressure measured here is the sum of the pre-column pressure, the pressure drop over the medium bed, and the post-column pressure. It is always higher than the pressure drop over the bed alone.

We recommend keeping the pressure drop over the bed below 1.5 bar. Setting the upper limit of your pressure gauge to 1.5 bar will ensure the pump shuts down before the medium is overpressured.

If necessary, post-column pressure of up to 3.5 bar can be added to the limit without exceeding the column hardware limit. To determine post-column pressure, proceed as follows:

Note: To avoid breaking the column, the post-column pressure must never exceed 3.5 bar.

Step	Action
1	Connect a piece of tubing in place of the column.
2	Run the pump at the maximum flow you intend to use for chromatography. Use a buffer with the same viscosity as you intend to use for chromatography. Note the back pressure as total pressure.
3	Disconnect the tubing and run at the same flow rate used in step 2. Note this back pressure as pre-column pressure.
4	Calculate the post-column pressure as total pressure minus pre-column pressure.

If the post-column pressure is higher than 3.5 bar, take steps to reduce it (shorten tubing, clear clogged tubing, or change flow restrictors) and perform steps 1–4 again until the post-column pressure is below 3.5 bar. When the post-column pressure is satisfactory, add the post-column pressure to 1.5 bar and set this as the upper pressure limit on the chromatography system.

First time use


Ensure an appropriate pressure limit has been set. Equilibrate the column for first time use or after long storage by running:

100 ml binding buffer e.g. PBS, pH 7.3 (140 mM NaCl, 2.7 mM KCl, 10 mM Na₂HPO₄, 1.8 mM KH₂PO₄, pH 7.3 at 5 ml/min.

GSTPrep FF 16/10 column can be used directly on ÄKTA design systems without the need for any extra connectors.

Try these conditions first

Binding buffer	PBS, pH 7.3 (140 mM NaCl, 2.7 mM KCl, 10 mM Na ₂ HPO ₄ , 1.8 mM KH ₂ PO ₄ , pH 7.3)
Elution buffer	50 mM Tris-HCl, 10 mM reduced glutathione, pH 8.0
Flow rate	Sample loading: 1–5 ml/min (30–150 cm/h) Washing and elution: 2–10 ml/min (60–300 cm/h)



Buffers and solvent resistance

De-gas and filter all solutions through a 0.45 µm filter to increase column lifetime.



Daily use

All commonly used aqueous buffers.



Cleaning

Guanidine hydrochloride, up to 6 M for 1 h at room temperature.
1 M acetate pH 4.0 for 1 h at room temperature
Ethanol, up to 70%.
Non-ionic detergents.




Avoid

Unfiltered solutions.


Sample preparations

Filter the sample through 0.45 µm filter or centrifuge at 10 000 × g for 10 min. If possible dissolve the sample in binding buffer.



Delivery/storage

The column is supplied in 20% ethanol. If the column is to be stored for more than two days after use, clean the column according to the procedure described under [Cleaning-in-place \(CIP\)](#), on page 2. Then equilibrate with at least 100 ml 20% ethanol at a flow rate of 5 ml/min.



Note: HiPrep columns cannot be opened or refilled.

Optimization

Perform your first run according to [Try these conditions first](#), on page 1. If the results are unsatisfactory, consider the following:

Action	Effect
Decrease flow rate	Increased binding capacity. Due to the relatively slow binding kinetics between GST and glutathione, it is important to keep the flow rate low during sample loading/elution for maximum binding capacity/elution.

The concentration of GST-tagged protein can be estimated by measuring the absorbance at 280 nm. The GST-tag can be approximated using the conversion; $A_{280} \approx 1$, which corresponds to ~0.5 mg/ml.

The concentration of GST-tagged protein may also be determined by standard chromogenic methods (for example Lowry, BCA, and Bradford assays). If Lowry or BCA assays are to be used, the sample must first be buffer exchanged using a HiTrap™ Desalting column, HiPrep 26/10 Desalting column, or dialyzed against PBS to remove glutathione, which can interfere with the protein measurement. The Bradford method can be used in the presence of glutathione.

The reuse of GSTPrep FF 16/10 depends on the nature of the sample and should only be performed with identical samples.

Cleaning-in-place (CIP)

If the medium appears to be losing binding capacity, it may be due to an accumulation of precipitate, denatured or non-specifically bound proteins.

Removal of precipitated or denatured substances:

- Wash with 40 ml of 6 M guanidine hydrochloride, immediately followed by 100 ml of PBS, pH 7.3 at a flow rate of 5 ml/min.

Removal of hydrophobically bound substances:

- Wash with 60–80 ml of 70% ethanol or 40 ml of detergent immediately followed by 100 ml of PBS, pH 7.3 at a flow rate of 5 ml/min.

Note: HiPrep columns cannot be opened or refilled.

Troubleshooting

Symptom	Remedy
Low yield of eluted protein from the column	<p>The 10 mM recommended in this protocol should be sufficient for most applications, but exceptions may occur. Try 50 mM Tris-HCl, 20–40 mM glutathione, pH 8.0 as elution buffer to increase efficiency of eluting the protein.</p> <p>A low pH may limit elution from GSTPrep FF 16/10. Increasing the pH of the elution buffer to pH 8–9 may improve elution without requiring an increase in the concentration of glutathione.</p> <p>Including 0.1–0.2 M NaCl to the elution buffer may improve results.</p> <p>Non-specific hydrophobic interactions may prevent solubilization and elution of tagged proteins from GSTPrep FF 16/10. Adding a non-ionic detergent may improve results.</p> <p>Adding 2% N-octylglucoside can significantly improve elution of some GST-tagged proteins.</p>
Increased back pressure	<p>Reverse the flow direction and pump 100 ml elution buffer through the column at a flow rate of 5 ml/min at room temperature. Return to normal flow direction and run 100 ml binding buffer at a flow rate of 5 ml/min. (Try different cleaning procedures described in Cleaning-in-place (CIP), on page 2.)</p>
Loss of resolution and/or decreased sample recovery	<p>Try different cleaning procedures described in Cleaning-in-place (CIP), on page 2.</p>
Air in the column	<p>Reverse the flow direction and pump 100 ml of well de-gassed binding buffer through the column at a flow rate of 5 ml/min at room temperature.</p>

Intended use

The GSTPrep FF 16/10 is intended for research use only, and shall not be used in any clinical or in vitro procedures for diagnostic purposes.

Ordering information

Product	No. per pack	Code No.
GSTPrep FF 16/10	1 x 20 ml	28936550

Related products		
GSTrap™ FF	2 x 1 ml	17513002

Related products		
GSTrap FF	5 x 1 ml	17513001
GSTrap FF	100 x 1 ml ¹	17513005
GSTrap FF	1 x 5 ml	17513101
GSTrap FF	5 x 5 ml	17513102
GSTrap FF	100 x 5 ml ¹	17513105
Glutathione	25 ml	17513201
Sepharose 4 Fast Flow		
Glutathione	100 ml	17513202
Sepharose 4 Fast Flow		
Glutathione	500 ml	17513203
Sepharose 4 Fast Flow		
HiPrep 26/10 Desalting	1 x 53 ml	17508701
HiPrep 26/10 Desalting	4 x 53 ml	17508702
HiTrap Desalting	5 x 5 ml	17140801
HiTrap Desalting	100 x 5 ml ¹	11000329

¹ Special pack size delivered on specific customer order.

Accessories		
HiTrap/HiPrep 1/16" male connector to ÄKTAdesign	8	28401081
To connect columns with 1/16" connections to FPLC System:		
Union M6 female/ 1/16" male	5	18385801

Related printed literature		
GST Gene Fusion System Handbook		18115758
The Recombinant Protein Purification Handbook, Principles and Methods		18114275
Affinity Chromatography Handbook, Principles and Methods		18102229
Affinity Chromatography Columns and Media, Selection guide		18112186
Prepacked chromatography columns for ÄKTAdesign and Ettan™ LC systems, Selection guide		28931778
Glutathione Sepharose, Selection guide		28916833

Further information

For more information, please refer to cytiva.com/protein-purification cytiva.com/purification_techsupport or refer to "GST Gene Fusion System Handbook" and "The Recombinant Protein Purification Handbook" which can be ordered, see Ordering information.

[cytiva.com/protein-purification](https://www.cytiva.com/protein-purification)

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.

Ettan, GSTprep, GSTrap, HiPrep, HiTrap and Sepharose are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

All other third-party trademarks are the property of their respective owners.

© 2020 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

28941327 AC V:4 10/2020

