

# SFM4HEK293

## HYCLONE MEDIA AND SUPPLEMENTS

HyClone™ SFM4HEK293 is a protein-free and animal-derived component-free (ADCF) cell culture medium. The product was developed through the HyClone Metabolic Pathway Design process (see box) to support the growth of human embryonic kidney (HEK) 293 cells and the production of adenoviral vectors and recombinant proteins. The medium is specifically formulated using our proprietary lipid complexing process for enhanced stability and broad-spectrum cell culture applications. The product is successfully tested in a variety of culture systems including T-flasks, shakers, roller bottles, and bioreactors. SFM4HEK293 is available in liquid and powder formats in user-friendly packaging (Fig 1).

Key features of SFM4HEK293 include

- Ready-to-use liquid formulation requires no supplementation
- Allows for direct or sequential adaptation
- Designed for large-scale culture applications, including perfusion and fed-batch strategies
- Manufactured from traceable components and according to cGMP guidelines

## Specifications

- Protein-free and ADCF formulation
- Liquid medium contains poloxamer 188
- Does not contain phenol red
- Developed through Metabolic Pathway Design process for high adenoviral vector and recombinant protein production

### Product handling

Store medium at 2°C to 8°C away from light. In addition, powder medium should be stored protected from moisture in a tightly sealed container.



**Fig 1.** SFM4HEK293 medium is available as liquid or powder in pack sizes suitable for small-volume cell culture as well as large-scale bioprocessing applications.

### Metabolic Pathway Design process

An optimal cell culture process is dependent on a variety of factors including the parental cell line, the genetic makeup of the specific clone, medium and feed composition, as well as process variables to maximize viable cell densities and titers while maintaining cell morphology. Our experts in medium design and development know and understand how these factors can influence the metabolic processes involved. They evaluate the culture's metabolic activities, measuring nutritional demand and waste creation to make sure the correct type and quantity of nutrients are used to minimize waste and resultant cell toxicity. Our experts use their understanding of metabolic pathways to optimize medium composition for enhanced productivity and viable cell densities. Once a medium has been optimized using this Metabolic Pathway Design process, our scientists can help you devise the most effective cell culture strategy using a combination of medium and feeds to further enrich productivity and reduce process inefficiencies.

## Suggested preparation

### Reconstitution of SFM4HEK293 powder medium

1. While stirring, add 18.4 g/L of SFM4HEK293 powder medium to cell culture-grade water (room temperature) at 90% of final preparation volume. If your water source is normally cool, it may be useful to adjust the water temperature. Using warmer room temperature water (22°C to 25°C) will improve solubilization time. Mix for 20 min or until dissolved.
2. Add 1.0 g/L poloxamer 188 and 2.0 g/L sodium bicarbonate. Minimum mixing time: 20 min.
3. Bring vessel to final volume with cell culture-grade water. Allow solution to mix for 10 to 20 min.
4. Check pH and osmolality. Expected values:
  - pH 7.0 to 7.4
  - Osmolality 290 to 340 mOsm/kg
5. Sterile filter into desired container using a 0.2 µm sterile filter.

### Preparation notes

SFM4HEK293 powder medium does not contain L-glutamine. Recommended concentration: 4 to 6 mM.

## General culture recommendations

1. Cultures should be incubated at 37°C in a 5% CO<sub>2</sub> environment.
2. The caps on culture flasks should be loosened and adequate vessel headspace should be given to provide gas exchange.
3. Seeding densities should be ~ 3.0 × 10<sup>5</sup> cells/mL. Cells should typically be subcultured every 3 to 5 days, as necessary.

### Direct adaptation

1. Transfer cells grown in current medium directly into SFM4HEK293 medium at 5.0 × 10<sup>5</sup> cells/mL.
2. When viable cell density reaches 2.0 to 4.0 × 10<sup>6</sup> cells/mL, subculture the cells.
3. Cells should be subcultured every 48 to 96 h.
4. If cell viability drops below 80%, proceed to sequential adaptation.

### Sequential adaptation

1. Transfer cells grown in current medium into SFM4HEK293 medium at a ratio of 1:1 using a seeding density of 5.0 × 10<sup>5</sup> cells/mL.
2. Incubate culture until two population doublings are observed. Subculture cells by mixing equal volumes of cell suspension in conditioned medium and fresh SFM4HEK293 medium (1:1 ratio).
3. Continue to subculture the cells using this method until the previously used medium is reduced below 0.05% concentration and cell viability is > 85%.

### Cryopreservation

SFM4HEK293 medium adapted cells can be cryopreserved in a 1:1 mixture of fresh and conditioned SFM4HEK293 medium supplemented with 10.0% DMSO.

## Quality control testing

Quality control test specifications are listed in Table 1.

**Table 1.** Test specifications<sup>1</sup>

|             |                               |
|-------------|-------------------------------|
| Appearance  | Clear solution                |
| Osmolality  | 290 to 340 mOsm/kg            |
| pH          | 7.0 to 7.4                    |
| Sterility   | No growth (bacteria or fungi) |
| Endotoxin   | < 10.0 EU/mL <sup>1</sup>     |
| Application | Growth promotion              |

<sup>1</sup> Refer to certificate of analysis for actual results.

## Custom production

Formulations and delivery systems can be customized to your specific process requirements or optimized to maximize process yields.

### Rapid Response Production (RRP)

Our Rapid Response Production (RRP) program manufactures up to 200 L of your custom prototype formulation within seven working days of your request. Use our RRP service to expedite the development and testing of custom buffers and process liquids for your biopharmaceutical manufacturing process.

**Table 2.** Supplement matrix

|                                 | Amino acids | Vitamins | Glucose | Trace elements | Growth factors | Hypoxanthine/thymidine | ADCF* lipids | ADCF cholesterol | Suitable for                           | Product code |
|---------------------------------|-------------|----------|---------|----------------|----------------|------------------------|--------------|------------------|--|--------------|
| Cell Boost 1 Supplement (R05.2) | •           | •        | •       |                |                |                        |              |                  | HEK293 CHO                             | SH30584      |
| Cell Boost 2 Supplement (R15.4) | •           |          | •       |                |                |                        |              |                  | PER.C6® CHO                            | SH30596      |
| Cell Boost 3 Supplement (JM3.5) | •           | •        | •       | •              |                | •                      |              |                  | Hybridoma Myeloma                      | SH30825      |
| Cell Boost 4 Supplement (PS307) | •           | •        | •       | •              | •              |                        | •            | •                | CHO                                    | SH30857      |
| Cell Boost 5 Supplement (CN-F)  | •           | •        | •       | •              | •              | •                      | •            | •                | Hybridoma NS0<br>HEK293 CHO            | SH30865      |
| Cell Boost 6 Supplement (CN-T)  | •           | •        | •       | •              | •              | •                      | •            | •                | T-Cells<br>Hybridoma NS0<br>HEK293 CHO | SH30866      |
| LS250 supplement                |             |          |         |                |                |                        | •            | •                | NS0                                    | SH30554      |
| LS1000 supplement               |             |          |         |                |                |                        |              | •                | NS0                                    | SH30555      |

\* Animal-derived component-free

## Related products

Table 2 gives an overview of HyClone supplements.

### HyClone Cell Boost kit

Cell Boost™ Process Supplements (100 g each) contain samples of supplements designed to increase cell productivity in a variety of cell lines (Table 2). Each supplement is developed through the Metabolic Pathway Design process and is chemically defined and protein-free with no animal derived components.

### HyClone LS250 supplement

LS250 is a chemically defined, animal-derived component-free lipid supplement developed to stimulate cell growth and MAb production in NS0 cell cultures using traditional hybridoma serum-free media.

### HyClone LS1000 supplement

LS1000 supplement is a chemically defined, animal-derived component-free lipid supplement developed to stimulate cell growth and MAb production in NS0 cell cultures using traditional hybridoma serum-free media. The supplement is formulated using a proprietary complexing process for enhanced cholesterol delivery. LS1000 has been successfully tested in a variety of serum-free medium cultures, including HyClone CDM4NS0 and CDM4MAb media.

## Ordering information

SFM4HEK293 is manufactured in homogenous liquid lot sizes up to 10 000 L and powder lots up to 250 000 L.

| Product   | Size           | Product code |
|---|----------------|--------------|
| HyClone SFM4HEK293 liquid medium<br>With L-glutamine    | 500 mL bottle  | SH30521.01   |
|   | 1000 mL bottle | SH30521.02   |
|   | 5 L bag        | SH30521.03   |
|   | 10 L bag       | SH30521.04   |
|   | 20 L bag       | SH30521.05   |
|   | 50 L bag       | SH30521.06   |
|   | 100 L bag      | SH30521.07   |
|   | 200 L bag      | SH30521.08   |
|   | 200 L bag      | SH30521.09   |
| HyClone SFM4HEK293 powder medium<br>Without L-glutamine | 1 × 5 L*       | SH30522.01   |
|   | 1 × 10 L*      | SH30522.02   |
|   | 1 × 50 L*      | SH30522.03   |
|   | 1 × 100 L*     | SH30522.04   |
|   | 1 × 500 L†     | SH30522.05   |
|   | 1 × 1000 L†    | SH30522.06   |

\* High-density polyethylene (HDPE) bottle

† Polybag/pail

| Related products                      | Size           | Product code |
|---------------------------------------|----------------|--------------|
| HyClone Cell Boost kit                | 6 × 100 g      | SH30890      |
| HyClone LS1000 cholesterol supplement | 50 mL bottle   | SH30554.01   |
|                                       | 100 mL bottle  | SH30554.02   |
|                                       | 500 mL bottle  | SH30554.03   |
|                                       | 1000 mL bottle | SH30554.04   |
| HyClone LS250 lipid supplement        | 100 mL bottle  | SH30555.01   |
|                                       | 500 mL bottle  | SH30555.02   |
|                                       | 1000 mL bottle | SH30555.03   |

## **[cytiva.com/hyclone](https://cytiva.com/hyclone)**

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. Cell Boost and HyClone are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

PER.C6 is a registered trademark of Crucell. All other third-party trademarks are the property of their respective owners.

© 2020 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit [cytiva.com/contact](https://cytiva.com/contact)

CY12974-10Jul20-DF

