

Hydrodynamic addressing in Biacore 4000 when using pre-immobilized sensor chips

This protocol describes how the hydrodynamic addressing (HA) procedure in Biacore 4000 should be run when using pre-immobilized sensor chips like Sensor Chip Protein A. For other sensor chips, the normal HA protocol should be run.

Preparation of solutions

1. Prepare running buffer HBS-EP+:


Mix 100 ml 10x HBS-EP+ (BR-1006-69) with 900 ml deionized water.

2. Prepare 10% HBS-EP+:

Mix 100 ml running buffer HBS-EP+ with 900 ml deionized water.

Hydrodynamic addressing

1. Dock a chip and select **Hydrodynamic Addressing with normalize**.
2. Place bottles with 10% HBS-EP+, deionized water and running buffer in the instrument and insert the tubing according to the picture below.



The only difference compared to the instructions in the software is that the tubing marked "Water" should be placed in the bottle with 10% HBS-EP+ solution:

10% HBS-EP+ solution in this position

Deionized water (as usual)

Running buffer HBS-EP+ (as usual)

Tray #1 Bottles

Position	Content	Estimated consumption
Water	Water	<100 ml, excluding bottle dead volume and standby consumption (450 ml/day)
A	Water	<100 ml, excluding bottle dead volume and standby consumption (110 ml/day)
B	Water	<100 ml, excluding bottle dead volume and standby consumption (10 ml/day)
1	Buffer	<100 ml, excluding bottle dead volume and standby consumption (75 ml/day)
2	Buffer	<100 ml, excluding bottle dead volume and standby consumption (75 ml/day)
3	Buffer	<100 ml, excluding bottle dead volume and standby consumption (75 ml/day)
4	Buffer	<100 ml, excluding bottle dead volume and standby consumption (75 ml/day)

3. Start the Hydrodynamic Addressing run.
4. **Important!** To avoid salt precipitation in the instrument, the 10% HBS-EP+ solution must be replaced with deionized water after the hydrodynamic addressing run and the tool **Change Buffers 1-4** should be run.



For local office contact information, visit
www.gelifesciences.com/contact

GE Healthcare UK Limited
Amersham Place
Little Chalfont
Buckinghamshire, HP7 9NA
United Kingdom

www.gelifesciences.com/biacore



GE, GE monogram and Biacore are trademarks of General Electric Company.

© 2016 General Electric Company
First published Jun. 2016

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare Bio-Sciences AB
Björkgatan 30, 751 84 Uppsala, Sweden

GE Healthcare Europe GmbH
Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare Bio-Sciences Corp.
100 Results Way, Marlborough, MA 01752, USA

GE Healthcare Dharmacon, Inc.
2650 Crescent Dr., Lafayette, CO 80026, USA

HyClone Laboratories, Inc.
925 W 1800 S, Logan, UT 84321, USA

GE Healthcare Japan Corporation
Sanken Bldg. 3-25-1, Hyakunincho Shinjuku-ku, Tokyo 169-0073, Japan