

# Creating a robust cell line development platform starts with the host cell line

Michelle Sabourin<sup>1</sup>, Paula Ravnikar<sup>1</sup>, Zhou Jiang<sup>1,2</sup>, Nan Lin<sup>2</sup>, and Peggy Lio<sup>1,2</sup>  
<sup>1</sup>Cytiva, 170 Locke Drive, Marlborough, MA 01752, USA; <sup>2</sup>Cytiva, 2780 N 200 W, North Logan, UT 84341, USA


**Abstract**

What do you need to create a robust cell line development platform? A robust host cell line – one that grows in suspension culture in the desired medium and has a consistent doubling time of < 24 h. We established such a host cell line, which can also seamlessly switch between cloning and production media.

Cytiva has adapted a serum-dependent CHO-K1 cell line to serum-free suspension culture in chemically defined and animal component-free HyClone™ media. To ensure success in the shortest time, we used a two-pronged, four-media adaptation strategy with eight adaptation trains. After growing the resulting GOCHO™ host cell line for 130 more generations, it doubled in < 20 h, displayed no cell clumping, and readily transitioned between cloning medium and HyClone™ ActiPro™ production medium. We then evaluated the GOCHO™ host cell line using model monoclonal antibodies. The cell line readily yielded stable clones that produced 3 to 5 g/L in a shake-flask, fed-batch process that was not optimized.

**Materials and methods**

CHO-K1 cells were thawed in α-MEM (Cytiva) + 10% fetal bovine serum (FBS; Cytiva). After a brief expansion in T-flasks, cells were seeded in shake flasks (SF) in media mixtures containing 50% of the desired end serum-free medium (SFM; Fig 1). Once success criteria were met at each stage, the cells were transitioned to the next media mixture containing less serum-containing medium (SCM). Cells were cryopreserved at each transition. Once cells were passaged into 100% SFM, a research cell bank (RCB) was generated and tested for mycoplasma. Cells were then thawed and passaged for an additional 100+ generations, with RCBs being generated every ~ 30 generations. Adaptation to multiple HyClone™ media was performed in parallel, with the goal of establishing a cell line that could transition seamlessly between cloning medium and the ultimate production medium, HyClone™ ActiPro™ (Cytiva). A seamless transition, defined as maintaining viability > 95% with no apparent growth lag, was tested first in SF (Fig 5), and later within the context of the cell line development workflow (see Fig 6 and 7 for representative clones).




**Summary**

CHO-K1 cells were adapted from adherent culture in SCM (Fig 2A) to suspension culture in serum-free chemically defined (Fig 2B) without clumping.


(A) (B)

**Fig 2.** Images of CHO-K1 (A) pre-adaptation and (B) post-adaptation.





**Seamless media switching**

One of the biggest challenges in cell line development (CLD) is the transition in/out of cloning media. To confirm that GOCHO™ cells were ready for CLD workflow testing, cells were grown in shake flasks and alternately passaged in either cloning media (CM) or ActiPro™ (AP) production media. The original train (green) was maintained in parallel for reference. Media switching was observed to be seamless as viability was maintained above 95% and there was no apparent growth lag. The differences in Td observed between CM and AP is consistent with the Td of GOCHO™ cells in CM and AP when not switching media.



**Fig 5.** Seamless media switching: GOCHO™ cells were thawed and alternately passaged between cloning medium and ActiPro™ medium. Doubling times (Td, bottom lines) and viability (circles) were monitored over 2 mo. Switches: 1 to CM, 2 to AP, 3 back to CM, and 4, back to AP.



**Conclusion**

CHO-K1 cells were successfully adapted out of adherent culture in SCM to suspension culture in SFM. The GOCHO™ host cell line was free of cell clumps (Fig 2), had a Td of < 20 h, and seamlessly switched between cloning media and ActiPro™ production media (Fig 5). GOCHO™ cell line performance within the context of the cell line development workflow was confirmed with the successful isolation of stable high-producing clones from three model monoclonal antibodies (Fig 6 and 7).